SORIMEX

PL Mankiety NIBP w.u.

EN Reusable NIBP cuffs

SORIMEX sp. z o.o. (dawniej SORIMEX sp. z o.o. sp. k.)
Równinna 25, 87-100 Toruń, Poland
office@sorimex.eu ♦ www.sorimex.eu

Rev. 2025-01-07

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PL

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EN

This instruction for use concerns:

- **Reusable NIBP cuffs** of the following type:
 - MSWZ - single-tube NIBP cuff with a bladder
 - MSWB - single-tube NIBP cuff bladderless
 - MDWZ - double-tube NIBP cuff with a bladder
 - MDWB - double-tube NIBP cuff bladderless

1. Intended use and operation of the product

NIBP cuffs are intended to be used for non-invasive blood pressure measurements by creating pressure on a relevant part of the patient's body.

The cuff is equipped with flexible connections, through which air is supplied via adapters, to create the required pressure.

NIBP cuffs are intended to be connected to a device monitoring patient's main vital functions.

2. Target patient groups

NIBP cuffs are intended to be used by qualified medical personnel in neonates, infants, children and adults.

3. Contraindications

No known contraindications to the use of the product.

4. Possible adverse effects

NIBP cuffs are devices having a direct contact with the patient's skin. There is a risk of irritation and allergy reactions in patients with particularly sensitive skin. If any skin reactions develop, medical advice should be sought.

5. Warnings

- 5.1. A damaged NIBP cuff should be disconnected and replaced with one in a good operating condition. NIBP cuffs must be used with units that are in a perfect operating condition.

6. Safety precautions

- 6.1. Do not use NIBP cuffs that are mechanically damaged.
- 6.2. Do not use NIBP cuffs that were in contact with any liquid.
- 6.3. Do not install the NIBP cuff on a limb on which the SpO2 sensor is installed.
- 6.4. Clean and disinfect the cuff before use.
- 6.5. Do not place the cuff on injured skin.
- 6.6. Do not leave the cuff on the limb after the monitoring has ended.
- 6.7. Excessive movements of a patient may result in measurements deviations.
- 6.8. Check a location for placing the cuff and the limb for signs of the obstructed blood flow, especially during a long-term monitoring.

7. How to use the cuff

- 7.1. Adjust the cuff to the circumference of the patient's arm/thigh
- 7.2. Install the cuff in such way that the  symbol is above the brachial/femoral artery.
- 7.3. Make sure that the cuff adheres to the patient's body, without excessive pressure and free space.
- 7.4. After the NIBP cuff is installed, ensure that the cuff hose is placed in a way reducing a possibility of patient's entangling or of damaging the device.

8. Cleaning and disinfection

These actions should be performed each time before the NIBP cuff is used. Before cleaning, remove the inner bladder from the cuff if it was installed in the cuff.

8.1. Cleaning

A cuff should be cleaned by delicately wiping it with a soft cloth wetted with warm water and soap or mild detergent. Before the use, the cuff should be dry.

During cleaning, do not wet or immerse the cuff, do not wash under running water, and do not use excessive force that may damage the cuff.

Do not use any corrosive agents or other aggressive chemical formulations. Avoid any contact with strong solvents: aromatic, chlorine, ketone, ethers or esters.

8.2. Disinfection

Disinfect the cuff by wiping with formulations intended for disinfection of medical devices made of plastic, if they do not negatively affect the cuff. Be sure to read the instructions for use of those formulations carefully before their use.

Materials recommended for disinfection of the cuff include 70% of isopropyl alcohol or 70% ethanol.

9. Signs of product degradation

Before each use, perform visual inspection of the product. NIBP cuffs cannot be used any longer when any of the below is noticed:

- integrity of the cuff surface is compromised;
- leaking of the cuff and/or its bladder;
- breaking of cuff air hoses;
- obstruction of cuff air hoses;
- hoses partially or completely torn out from the cuff

Any of the above defects disqualifies the product from further use. A product with such defects should either be repaired or disposed of.

10. Storage conditions

10.1. Storage conditions:

- temperature: 0°C to + 40°C,
- protect against sunlight, moisture and contamination

10.2. NIBP cuffs should be stored in the specified storage conditions, to prolong the product life as far as possible.

10.3. During storage, NIBP cuffs should be loosely wound, without any folds on the surface.

11. Compatibility

The NIBP cuff operates in connection with a device via a relevant NIBP adapter (connector).

NIBP adapters are components connecting NIBP cuffs to devices used for non-invasive blood pressure measurements.

Adapters include connectors and air hoses.

A type of an air hose and connectors should match a relevant design of the device with which it is to work.

12. Warranty

The product life is 36 years of the production date specified in the label. The guarantee period is 12 months of the sales date.

The general guarantee terms and conditions are available at the manufacturer website, www.sorimex.pl. The guarantee does not cover defects resulting from a failure to observe provisions of this operating manual.

13. Description of symbols used

Catalogue number

Pieces per packaging

Batch code

Medical device

Date of manufacture

Manufacturer

Consult instructions for use or consult electronic instructions for use

Caution

Keep away from sunlight

CE mark (conformance to the Regulation EU 2017/745)

Keep dry

Temperature limit

Non-sterile

Do not contain latex

Humidity limitation

14. Notification of incidents related to the product use

Each serious incident related to the product should be notified to the manufacturer and to an authority of the Member State competent for the user or patient's place of residence.

15. Product disposal

The used cuff should be disposed of in accordance with internal hygiene procedures in force at a given health care facility and with local regulations.