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Shenzhen Sino-K Medical Technology Co., Ltd
Room 401, Bldg2, Veteran Ind City, Gonghe Community
Xixiang Street, Baoan District, Shenzhen, Guangdong, China
Tel: 86-755-23341807, Fax: 86-755-23204076
info@sino-k.com

Instruction for Use

1. Intended use: The spo2 sensor is connected with the compatible pulse oximeter for continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate(PR) for adult patients weighing greater than 40 kg, pediatric patients weighing 10-50 kg, infant weighing 3-15 kg and neonatal weighing 1-3 kg. The SpO2 sensor is intended for use by qualified medical personnel.

2. Models

Usage	Category	Intended person	Measurement part	Proximal connector	Compatibility
Reusable	SP9 Serial	Adult (>40 kg)	Fingers or toes	Finger clip	See the label
	SP8 Serial	Adult (>40 kg)	Fingers or toes	Soft tip	See the label
	SP7 Serial	Pediatric (10-50 kg)	Fingers or toes	Finger clip	See the label
	SP6 Serial	Pediatric (10-50 kg)	Fingers or toes	Soft tip	See the label
	SP5 Serial	Neonate (<3 kg)	Palm or sole	Wrap	See the label
	SP4 Serial	All	Palm or sole	Multi-site Y	See the label
	SP3 Serial	Adult (>40 kg)	Ear	Ear clip	See the label
	SP18 Serial	Neonatal (<3 kg)/Adult (>40 kg)	Palm or sole / Fingers or toes	Sponge	See the label
	SP19 Serial	Pediatric (10-50 kg)	Fingers or toes	Sponge	See the label
	SP10 Serial	Neonatal (<3 kg)/Adult (>40 kg)	Palm or sole / Fingers or toes	Nonwovens	See the label
Disposable	SP11 Serial	Infant (3-20 kg)	Palm or sole	Nonwovens	See the label
	SP12 Serial	Pediatric (10-50 kg)	Fingers or toes	Nonwovens	See the label
	SP13 Serial	Adult (>40 kg)	Fingers or toes	Nonwovens	See the label
	SP14 Serial	Neonatal (<3 kg)/Adult (>40 kg)	Palm or sole / Fingers or toes	Foam	See the label
	SP15 Serial	Infant (3-20 kg)	Palm or sole	Foam	See the label
	SP16 Serial	Pediatric (10-50 kg)	Fingers or toes	Foam	See the label
	SP17 Serial	Adult (>40 kg)	Fingers or toes	Foam	See the label

Note: The SpO2 sensor is designed to match the specifications of the original equipment manufacturer, therefore ensure the appropriate sensor model number is selected for a particular pulse oximeter technology. Each sensor has unique labeling and specifications designed for compatibility with the specific monitor manufacturer (Nellcor, Nonin, BCI, GE and so on).
Note: See the compatibility labeling of SpO2 sensor, and choose the proper devices matching the sensor.

3. Technical description

SpO2 Accuracy (NO motion)	±2 % (70 - 100 %)
Pulse rate Accuracy (NO motion)	±3 (20 - 250 bpm)
Wavelengths and output Power	RED: 660 - 666 nm, IR: 940/905 nm, Output Power: less than 20mW

4. Working principle: The SpO2 measurement is based on the absorption of pulse blood oxygen to red and infrared light by means of sensor and SpO2 measuring unit(Pulse Oximeter). The light-electronic transducer in sensor converts the pulse red and infrared light modulated by pulse blood oxygen into electrical signal, the signal is processed by hardware and software of SpO2 measuring unit(Pulse Oximeter). The pleth curve or (and) numeral value of SpO2 will be obtained.

5. Contraindications: The spo2 sensor is contraindicated for use on active patients or for prolonged use.

6. Summary of SpO2 & pulse rate accuracy testing

- According to ISO80601-2-61, the 250 data pairs (SPO2 VS SaO2) were obtained throughout the hypoxia breathe-down trials in all 10 volunteers. The SpO2 accuracy is ±2% which meets the requirement of ISO80601-2-61.
- The pulse rate accuracy is less than ±3 in range 20~245 bpm, covering the entire pulse rate display range of 20~250bpm, which meets the requirements of ISO80601-2-61.

7. Instructions for use

- Hold the sensor with its opening towards the patient's index finger, toes, palm or sole. The side with a finger or toes sign is placed on the top.
- Put the patient's measurement part into the sensor until the tip against at the end of the sensor
- Adjust the measurement part to be placed evenly on the middle base of the sensor.
- Direct the cable along the top of the patient's hand or foot.
- Plug the sensor into the proper oximeter and verify operation as described in the oximeter operator's manual.
- Inspect the sensor site every 2 hours for skin integrity.

8. Operation and storage environment

Temperature: +5 to+40°C, Humidity: ≤80% (non-condensing), Atmospheric pressure: 86kPa-106kPa, Power supply: 3-5V DC
 Temperature: -10 to+40°C, Humidity: ≤80% (non-condensing), Atmospheric pressure: 86kPa-106kPa

9. WARNING: When disconnecting, pull on the plug itself, not on the cable

- This sensor is only use with compatible patient monitors or pulse oximeter.
- User should verify the compatibility of the monitor, sensor, and cable before use, otherwise cause patient injury.
- Check the site every 2 hours (more frequently if perfusion is poor).
- Routinely check to ensure adequate distal circulation to the sensor site.
- Carefully route cables to reduce any possibility of patient entanglement or strangulation.
- Patient conditions (such as reddening, blistering, skin discoloration, ischemic skin necrosis, and skin erosion) may warrant changing the site frequently or using a different style of sensor.
- Do not use the sensor if the sensor or the cable appears damaged.
- This device is not intended for use in a magnetic resonance imaging (MRI) environment.
- Conditions that may cause inaccurate reading and impact alarms include interfering substances, excessive ambient light, electrical interference, excessive motion, low perfusion, low signal strength, incorrect sensor placement, poor sensor fit, and movement of the sensor on the patient.
- Do not use a blood pressure cuff or arterial blood pressure measurement device on the same limb as the sensor site.
- It is possible for any device to malfunction; therefore, always verify unusual data by performing a formal patient assessment.
- The sensor should be disinfected before each use.
- Disposal of the sensors should comply with local regulation.

10. Cleaning & Disinfection

- Unplug the sensor before cleaning or disinfecting.
- Clean the part of the sensor that comes into contact with the skin by immersing it in a mild detergent solution or a 70% isopropyl alcohol solution. If low-level disinfection is required, use a 1:10 bleach solution. Take care to avoid contact of the plug with the cleaning agent.
- Clean the sensor part (not the connector) with the wet cloth, and then dry it up with the dry cloth in the end.

11. Reporting incidents involving the device: Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident.

12. Warranty: Offers a one-year warranty against defects in material and /or workmanship from the date of purchase.

13. Shelf Life: 3 years

14. Duration of use: recommendations for every 12 months to replace once. If it is necessary to use, please to evaluate, and the use of period should be tightened inspection.

15. Symbol: The following symbols may appear on the product or product labeling:

	Catalogue number		Unique Device Identifier		Protect from moisture		Do not reuse
	Quantity in package		Manufacturer		Labeling in package and electronic equipment		Humidity limit
	Batch number		Authorized representative in the EU		Warning		Atmospheric pressure limit
	Expire date		Importer		Protect from sunlight		Non-sterile
	Date of manufacture		Medical device		Product does not contain latex		Refer to instruction manual
	CE conformity marking		Do not use if package is damaged		Temperature limit		IPX2 Water protection

SORIMEX

SORIMEX sp. z o.o.
 Równina 25, 87-100 Toruń, Poland
 office@sorimex.eu ♦ www.sorimex.eu
 Prolinx GmbH
 Brehmstr. 56, 40239, Duesseldorf
 Tel: 49-2131-4051968-0, Fax: 49-2131-4051968-9
 info@eulinx.eu

Shenzhen Sino-K Medical Technology Co.,Ltd
 Room 401, Bldg2, Veteran Ind.City, Gongfe Community
 Xixiang Street, Baoan District, Shenzhen,
 Guangdong, China
 Tel: 86-755-23341807, Fax: 86-755-23204076
 info@sino-k.com

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